



Rabbit anti-Goat IgG Secondary
Antibody(unconjugated)

#L3041

Catalog Number:

L3041-1, L3041-2

Amount:

Bulk Order, 1mg (IgG content: 10mg/mL)

Form of Antibody:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl.

Storage/Stability:

Store at -20°C/1 year

Immunogen:

The antiserum was developed in Rabbit using the GoatIgG as the immunogen.

Purification:

The antibody is obtained by affinity purified technique, which removes essentially all Rabbit serum proteins, except the specific antibody for Goat IgG.

Specificity/Sensitivity:

The antibody conjugate is specific for Goat IgG. No cross-reactivity was observed with Goat IgA, IgM proteins by ELISA.

Identification:

The purity of IgG was $\geq 95\%$ basing on the SDS-PAGE analysis.

The protein content of IgG was **10mg/ml** tested by Lowry method.