



## Hexokinase 1

### Mouse monoclonal Antibody

#53543

**Catalog Number:** 53543

**Amount:** 100µg/100µl

**Swiss-Prot No. :** P19367

**Gene name:**hk1

**Gene id:**3098

**Clone Number:** 6B8-G11-E7

**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol

**Storage/Stability:** Store at -20°C/1 year


**Immunogen:** Purified recombinant human Hexokinase 1 protein fragments expressed in E.coli

**Purification:** affinity-chromatography

**Specificity/Sensitivity:** This antibody detects endogenous levels of Hexokinase 1 and does not cross-react with related proteins

**Reactivity:** Human, Mouse, Rat

**Applications:** Predicted MW: 102kd WB: 1:1000



Western blot detection of Hexokinase 1 in MCF7, Jurkat, CHO-K1, C2C12, C6, 3T3 and HeLa cell lysates using Hexokinase 1 mouse mAb (1:1000 diluted). Predicted band size: 102KDa. Observed band size: 102KDa.

### Background:

Hexokinases phosphorylate glucose to produce glucose-6-phosphate, the first step in most glucose metabolism pathways. This gene encodes a ubiquitous form of hexokinase which localizes to the outer membrane of mitochondria. Mutations in this gene have been associated with hemolytic anemia due to hexokinase deficiency. Alternative splicing of this gene results in five transcript variants which encode different isoforms, some of which are tissue-specific. Each isoform has a distinct N-terminus; the remainder of the protein is identical among all the isoforms. A sixth transcript variant has been described, but due to the presence of several stop codons, it is not thought to encode a protein.