



SMAD3 (1E8)


Mouse monoclonal Antibody

#55038

Catalog Number: 55038**Amount:** 100µg/100µl**Swiss-Prot No. :** P84022**GENE ID:** 4088**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of SMAD3 and does not cross-react with related proteins**AlternativeNames:** Mothers against decapentaplegic homolog 3 (MAD homolog 3) (Mad3) (Mothers against DPP homolog 3) (hMAD-3) (JV15-2) (SMAD family member 3) (SMAD 3) (Smad3) (hSMAD3)**Reactivity:** Human, Mouse, Rat**Applications:**

Predicted MW: 52kd

IHC: 1:100 -200



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma
Tissue using Smad3 Mouse mAb

Background: