



Collagen I (2D10)

Mouse monoclonal Antibody


#55203

Catalog Number: 55203**Amount:** 100µg/100µl**Swiss-Prot No. :** P02452**GENE ID:** 1277/1278**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of Collagen I and does not cross-react with related proteins**Alternative Names:** Collagen alpha-1(I) chain (Alpha-1 type I collagen)**Reactivity:** Human, Mouse, Rat

Applications:

Predicted MW: 139kd

IHC: 1:50-200



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma

Tissue using Collagen I Mouse mAb

Background: