

TBC1D3



Catalog Number: 24115-1, 24115-2

Amount: 50µg/50µl, 100µg/100µl

Swiss-Prot No. :Q8IZP1/A6NDS4/Q6IPX1

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl,0.02% sodium azide and 50% glycerol.

Storage/Stability: Store at -20°C/1 year

Immunogen: The antiserum was produced against synthesized peptide derived from Human TBC1D3 Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.


Specificity/Sensitivity: TBC1D3 antibody detects endogenous levels of total TBC1D3 protein

Reactivity: Human

Applications:

Predicted MW: 62kd

WB:1:500~1:1000



Western blot analysis of extracts from HeLa cells, using TBC1D3 antibody.

Background :