




PBOV1 Antibody

#24116

Catalog Number: 24116-1, 24116-2**Amount:** 50µg/50µl, 100µg/100µl**Swiss-Prot No. :** Q9GZY1**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.**Storage/Stability:** Store at -20°C/1 year**Immunogen:** The antiserum was produced against synthesized peptide derived from Human PBOV1**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.**Specificity/Sensitivity:** PBOV1 antibody detects endogenous levels of total PBOV1 protein**Reactivity:** Human**Applications:**

Predicted MW: 19kd

WB:1:500~1:1000



Western blot analysis of extracts from K562 cells, using PBOV1 antibody.

Background :