



CO4A2 Antibody

#24118


Catalog Number: 24118-1, 24118-2**Amount:** 50µg/50µl, 100µg/100µl**Swiss-Prot No. :** P08572**Form of Antibody:** Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.**Storage/Stability:** Store at -20°C/1 year**Immunogen:** The antiserum was produced against synthesized peptide derived from Human CO4A2**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.**Specificity/Sensitivity:** CO4A2 antibody detects endogenous levels of total CO4A2 protein**Reactivity:** Human, Mouse**Applications:**

Predicted MW: 150kd

WB: 1:500~1:1000

IHC 1:50:200

IF/ICC: 1:100-500



Western blot analysis of extracts from Jurkat cells, using Collagen IV α 2 antibody.

Background :