

## PMEPA1

Order: order@swbio.com

**Catalog Number:** 24117-1, 24117-2 **Amount:** 50µg/50µl, 100µg/100µl Swiss-Prot No.: Q969W9

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM

NaCl,0.02% sodium azide and 50% glycerol. Storage/Stability: Store at -20℃/1 year


Immunogen: The antiserum was produced against synthesized peptide derived from Human PMEPA1 Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity/Sensitivity: PMEPA1 antibody detects endogenous levels of total PMEPA1 protein

Reactivity: Human

**Applications:** 

Predicted MW: 37kd WB:1:500~1:1000 IHC1:50:200 IF/ICC:1:100-500



Western blot analysis of extracts from HT-29 cells, using TMEPA antibody.

## Background: