

## PIAS4

Order: order@swbio.com

Catalog Number: 24120-1, 24120-2 Amount: 50µg/50µl, 100µg/100µl Swiss-Prot No.: Q8N2W9

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM

NaCl,0.02% sodium azide and 50% glycerol. Storage/Stability: Store at -20°C/1 year


Immunogen: The antiserum was produced against synthesized peptide derived from Human PIAS4 Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Specificity/Sensitivity: PIAS4 antibody detects endogenous levels of total PIAS4 protein

Reactivity: Human, Mouse

**Applications:** 

Predicted MW: 57kd WB:1:500~1:2000 IHC1:50:200



Western blot analysis of PIAS4 expression in Jurkat cell extract

## Background:

Functions as an E3-type small ubiquitin-like modifier (SUMO) ligase, stabilizing the interaction between UBE2I and the substrate, and as a SUMO-tethering factor. Plays a crucial role as a transcriptional coregulation in various cellular pathways, including the STAT pathway, the p53 pathway, the Wnt pathway and the steroid hormone signaling pathway.