

VEGF-A

Order: order@swbio.com

Catalog Number: 25516-1, 25516-2 Amount: 50µg/50µl, 100µg/100µl Swiss-Prot No.: P15692

Form of Antibody: Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM

NaCl,0.02% sodium azide and 50% glycerol. Storage/Stability: Store at -20°C/1 year

Immunogen: .Synthesized peptide derived from the Internal region of human VEGF-A.

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using


epitope-specific immunogen.

Specificity/Sensitivity: VEGF-A Antibody detects endogenous levels of total VEGF-A

Reactivity: Human


Applications:

Predicted MW:93kd WB:1:500-2000 IHC:1:100-300



Western Blot analysis of 22RV-1, SW480, BT474, K562,

mouse liver cells using VEGF-A Polyclonal Antibody.



Immunohistochemical analysis of paraffin-embedded human-heart

Background: