



Goat anti-Chicken IgG Secondary
Antibody(FITC conjugated)

#L3205

Catalog Number:

L3205-1, L3205-2

Amount:

Bulk Order, 1ML

Form of Antibody:

FITC-labeled Goat IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl.

Storage/Stability:

Store at -20°C/1 year

Immunogen:

The antiserum was developed in Goat using the Chicken IgG as the immunogen.

Purification:

The antibody is obtained by affinity purified technique, which removes essentially all Goat serum proteins, except the specific antibody for Chicken IgG. Goat anti-Chicken IgG is conjugated to FITC by means of a method described by Wilson MB and Nakane PK.

Specificity/Sensitivity:

The antibody conjugate is specific for Chicken IgG. No cross-reactivity was observed with Chicken IgA, IgM proteins by ELISA.

Applications:

Flow-Cyt/ IF=1:250-2500

References:

Wilson MB and Nakane PK. In Immunofluorescence and Related Staining Techniques, Elsevier/North Holland Biomedical Press, Amsterdam, P215 (1978).