



Oct

Mouse monoclonal Antibody

#55250


**Catalog Number:** 55250**Amount:** 100µg/100µl**Swiss-Prot No. :** P14859**GENE ID:** 5451**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of Oct and does not cross-react with related proteins**Alternative Names:** XRCC4; DNA repair protein XRCC4; X-ray repair cross-complementing protein 4**Reactivity:** Human**Applications:**

Predicted MW: 89kd

WB:1:500-2000

IHC:1:200

IF:1:100-300



Western blot analysis of 1) Hela, 2) Jurkat, 3) HepG2,

**Background:**