



MAP2


Mouse monoclonal Antibody

#55243

Catalog Number: 55243**Amount:** 100µg/100µl**Swiss-Prot No. :** P11137**GENE ID:** 4133**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of MAP2 and does not cross-react with related proteins**Alternative Names:** Microtubule-associated protein 2; MAP-2**Reactivity:** Human, Mouse, Rat

Applications:

Predicted MW: 20kd IHC:1:200 IF:1:100-200



IHC staining of Human brain tissue paraffin-embedded