



Cytokeratin 6

Mouse monoclonal Antibody

#55246


Catalog Number: 55246**Amount:** 100µg/100µl**Swiss-Prot No. :** P02538/P04259/P48668**GENE ID:** 3854/286887/3853**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of Cytokeratin 6 and does not cross-react with related proteins**AlternativeNames:****Reactivity:** Human**Applications:**

Predicted MW: 56-60kd

WB:1:1000

IHC:1:200

IF:1:100-300



Western blot analysis of 1) Hela, 2) HepG2

Background: