



IDE


Mouse monoclonal Antibody

#55247

Catalog Number: 55247**Amount:** 100µg/100µl**Swiss-Prot No. :** P14735**GENE ID:** 3416**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of IDE and does not cross-react with related proteins**Alternative Names:** IDE; Insulin-degrading enzyme; Abeta-degrading protease; Insulin protease; Insulinase; Insulysin**Reactivity:** Human**Applications:**

Predicted MW: 118kd

WB:1:1000 IHC:1:200 IF:1:100-300



Western blot analysis of 1) Hela, 2) HepG2

Background: