


**XRCC4****Mouse monoclonal Antibody****#55249****Catalog Number:** 55249**Amount:** 100µg/100µl**Swiss-Prot No. :** Q13426**GENE ID:** 7518**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of XRCC4 and does not cross-react with related proteins**Alternative Names:** XRCC4; DNA repair protein XRCC4; X-ray repair cross-complementing protein 4**Reactivity:** Human**Applications:**

Predicted MW: 38kd WB:1:2000 IHC:1:200 IF:1:100-300 IP:1:200

**Background:**