



**Tau****Mouse monoclonal Antibody****#55260****Catalog Number:** 55260**Amount:** 100µg/100µl**Swiss-Prot No. :** P10636**GENE ID:** 4137**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of Tau and does not cross-react with related proteins**AlternativeNames:** MAPT**Reactivity:** Human, Mouse, Rat**Applications:**

Predicted MW: 79kd

WB: 1:1000-2000

IHC: 1:50-200



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue

Background: