


**FH Fumarase****Mouse monoclonal Antibody**

#55257

**Catalog Number:** 55257**Amount:** 100µg/100µl**Swiss-Prot No. :** P07954**GENE ID:** 2271**Form of Antibody:** Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:** This antibody detects endogenous levels of FH Fumarase and does not cross-react with related proteins**AlternativeNames:****Reactivity:** Human, Mouse, Rat**Applications:**

Predicted MW: 50kd      WB: 1:1000-2000      IF: 1:50-200



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain tissue,

**Background:**