




Ubiquitin

Mouse monoclonal Antibody

#55261

Catalog Number: 55261**Amount:** 100µg/100µl**Swiss-Prot No. :**PAN**GENE ID:**7314**Form of Antibody:**Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:**This antibody detects endogenous levels of Ubiquitin and does not cross-react with related proteins**AlternativeNames:****Reactivity:** Human, Mouse, Rat**Applications:**

Predicted MW: N/A WB:1:1000-2000 IHC:1:50-200



Western blot analysis of 1) HeLa Cell Lysate, 2) 3T3 Cell Lysate, 3) Rat Brain Tissue

Background: