


**AKT(P-Ser463)****Mouse monoclonal Antibody****#55263****Catalog Number:** 55263**Amount:** 100µg/100µl**Swiss-Prot No. :**P31749**GENE ID:**207**Form of Antibody:**Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol**Storage/Stability:** Store at -20°C/1 year**Purification:** affinity-chromatography**Specificity/Sensitivity:**This antibody detects endogenous levels of AKT(P-Ser463) and does not cross-react with related proteins**AlternativeNames:** AKT1**Reactivity:** Human**Applications:**

Predicted MW: 60      WB:1:1000-2000      IHC:1:50-200



Western blot analysis of PC3 Cell Lysate

**Background:**